## In the Claims:

1. (Currently Amended) A method of diagnosing or monitoring a lysosomal storage disorder in a patient subject, comprising:

obtaining a first sample from the patient subject;

measuring a first level of at least a first saposin in the first sample obtained from the patient subject;

comparing the first level to a baseline level, wherein the baseline level is the level of at least the first saposin as determined in a control population of <u>patients-subjects</u> unaffected by the lysosomal storage disorder; and

determining <u>an absence</u>, a presence or extent of a lysosomal storage disorder when the first level is similar or different than the 95<sup>th</sup> percentile of the baseline level of at least the first saposins in the control population; wherein,

- (i) the similarity of the first level compared to the baseline level is an indicator of absence of the lysosomal storage disorder in the patient-subject;
- (ii) the difference of the first level compared to the baseline level is an indicator of presence or extent of the lysosomal storage disorder in the patient-subject;

4

- (iii) the first saposin comprises saposin A, saposin B, saposin C, saposin D, or a combination thereof; and
- (iv) the first sample is a plasma, serum, <u>or</u> whole blood<del>, urine, or amniotic fluid</del> sample; <u>and</u>
- (v) the lysosomal storage disorder is selected from the group consisting of Cystinosis, Fabry disease, Galactosialidosis, Gaucher disease, GM I-gangliosidosis, I-cell disease, Krabbe disease, a-Mannosidosis, Metachromatic Leukodystrophy, MPS I, MPS II, MPS IIIA, MPS IIIB, MPS IIIC, MPS IIID, MPS IVA, MPS VI, Multiple Sulphatase Deficiency, Neuronal Ceroid Lipofuscinoses, Niemann-Pick disease (A/B), Niemann-Pick disease (C), Pompe disease, Sandhoff disease, Sialic Acid Storage disease, Tay-Sachs disease type I, Tay-Sachs disease (A/B), and Wolman disease.
- 2. Cancelled.

- 3. Cancelled.
- 4. (Previously Amended) The method of claim 1, further comprising indicating a presence of the lysosomal disorder in the patient-subject when the first level exceeds the baseline level.
- 5. (Currently Amended) The method of claim 1, further comprising:

measuring a second level of a second saposin in a second sample from the patient subject, wherein the first saposin and second saposin are the same, and the first and second samples are obtained at different times; and

comparing the first level and the second level in the samples to monitor progression of the disease,

determining <u>an absence</u>, a presence or extent of a lysosomal storage disorder when the second level is similar or different than the 95<sup>th</sup> percentile of the baseline level of at least the two saposins in the control population;

wherein,

- (i) the second saposin comprises saposin A, saposin B, saposin C, saposin D, or a combination thereof;
- (ii) the comparison of the first level and the second level is an indicator of the progression of the disease in the patient subject; and
- (iii) the second sample is a plasma, serum, <u>or</u> whole blood<del>, urine, or amniotic</del> fluid sample.
- 6. (Previously Amended) The method of claim 1, further comprising selecting the patient subject that is undergoing treatment for the lysosomal storage disorder.
- 7. Cancelled.
- 8. (Previously Amended) The method of claim 1, further comprising selecting the patient subject that is not known to have a lysosomal storage disorder before the measuring step.
- 9. (Previously Amended) The method of claim 1, further comprising selecting the patient subject that is an infant less than one year old.

5

- 10. (Previously Amended) The method of claim 1, further comprising selecting the patient subject that is a fetus and the sample is a fetal blood sample.
- 11. (Previously Amended) The method of claim 5, wherein a change in the first level of the saposin indicates progression or regression of the disorder in the patient subject that is known to have a lysosomal storage disorder.
- 12. (Previously Amended) The method of claim 5, wherein a change in the first level of the saposin indicates a response to treatment of the lysosomal storage disorder in the patient subject that being treated for the lysosomal storage disorder.
- 13. Cancelled.
- 14. Cancelled.
- 15. (Previously Presented) The method of claim 1, wherein the measuring step comprises detecting binding between a saposin polypeptide and an antibody.
- 16. (Original) The method of claim 15, wherein the antibody is a monoclonal antibody.
- 17. (Original) The method of claim 15, wherein the antibody is immobilized to a solid phase.
- 18. (Previously Amended) The method of claim 1, wherein the lysosomal storage disorder is selected from the group consisting of cystinosis, Fabry's disease, Niemann-Pick disease, Pompe's disease, Wolman disease.
- 19. (Previously Amended) The method of claim 1, further comprising informing the patient subject or a parent or guardian thereof of the presence of the lysosomal storage disorder.
- 20. (Previously Presented) The method of claim 1, further comprising determining a treatment program based on the measurement of the first level of the first saposin.
- 21. (Withdrawn) A method of diagnosing or monitoring a lysosomal storage disorder in a patient, comprising: measuring the level of a-glucosidase in a tissue sample from a patient, wherein the level is an indicator of the presence or extent of the disorder in the patient.
- 22. (Withdrawn) The method of claim 21, wherein the sample is a plasma sample.

ATTORNEY DOCKET: MAYO-0005 (128675.00014)

- 23. (Withdrawn) The method of claim 21, wherein the sample is a blood sample.
- 24. (Withdrawn) The method of claim 21, further comprising diagnosing the presence of a disorder selected from the group consisting of acid lipase disease, mannosidosis, MPSII, MPS IIIA, MSD, mucolipidosis, N-P (A/B), N-P (C), Sandhoff, SAS or TSD B1, if the measured level of a-glucosidase exceeds the mean level in a control population of individuals not having a lysosomal storage disease.
- 25. (Withdrawn) The method of claim 21, further comprising diagnosing the presence of disorder selected from the group consisting of galactosialidosis, MPS IVA and Pompe's disease if the measured level of a-glucosidase is below the mean level in a control population of individuals not having lysosomal storage disease.
- 26. (Withdrawn) A method of diagnosing a lysosomal storage disorder comprising measuring a level of a saposin in a tissue sample from the patient; measuring a level of LAMP-1 or LAMP-2 in a second tissue sample from the patient; measuring a level of a glucosidase in a third tissue sample from the patient; wherein an increased level of saposin and/or LAMP-1 or LAMP-2, and/or an increased or decreased level of a-giucosidase in the sample relative to respective mean levels in a control population is an indicator of presence or extent of the disorder in the patient.
- 27. (Withdrawn) A method of diagnosing Pompe's disease in a patient, comprising measuring a level of a saposin in a tissue sample from the patient; measuring the level of a-glucosidase in a second tissue sample from the patient; wherein the presence of an increased level of the saposin and a decreased level of the a-glucosidase relative to mean levels of the saposin and a-glucosidase in a control population of individuals not having a lysosomal storage disorder indicates Pompe's disease or susceptibility thereto.
- 28. (Withdrawn) A method of screening patients for presence of lysosomal storage disorder, comprising: measuring the level of a LAMP-1 polypeptide in a sample from the patient; measuring the level of a saposin peptide in the sample, the presence of an increased level of LAMP-1 or saposin or both relative to mean levels in a control population, indicating susceptibility to a lysosomal disorder.

7

- 29. (Withdrawn) A diagnostic kit comprising: a first reagent that binds to a LAMP; a second reagent that binds to a saposin.
- 30. (Withdrawn) The diagnostic kit of claim 29, further comprising a third reagent that binds to a glucosidase.
- 31. (Withdrawn) The diagnostic kit of claim 30, wherein the first, second and third reagents are antibodies.
- 32. (Withdrawn) In a method of screening a patient for presence or susceptibility to disease, comprising performing a plurality of diagnostic tests on a tissue sample from the patient for a plurality of diseases, the improvement wherein one of the diagnostic tests comprises measuring the level of a saposin.
- 33. (Withdrawn) In the method of claim 32, the further improvement wherein a second of the diagnostic tests comprising measuring the level of LAMP-1 in the tissue sample from the patient.
- 34. (Withdrawn) In the method of claim 33, the further improvement wherein a third of the diagnostic tests comprises measuring the level of a-glucosidase in the tissue sample from the patient.
- 35. (Withdrawn) In the method of claim 32, the further improvement wherein a fourth of the diagnostic test comprises analysing a nucleic acid encoding an enzyme associated with a lysosomal storage disorder for a polymorphic form correlated with the disorder.
- 36. (Previously Amended) A method of monitoring treatment of a lysosomal storage disease in a patient, comprising:

determining a pre-treatment baseline level of a saposin in a sample from the patient with a lysosomal storage disorder before treatment with an agent;

determining a post-treatment baseline level of the saposin in a sample from the patient with the lysosomal storage disorder after treatment with the agent; and

comparing the pre-treatment baseline level of the with the post-treatment baseline level of the saposin, wherein

8

- (i) the sample is a plasma, serum, whole blood, urine, amniotic fluid sample, or a mixture of;
- (ii) saposin is selected from the group consisting of saposin A, saposin B, saposin C, saposin D, and
- (iii) a reduction in the post-treatment baseline level relative to the pre-treatment baseline level indicates a positive treatment outcome.
- 37. (Withdrawn) A method of monitoring treatment of acid lipase disease, mannosidosis, MPSII, MPS IIIA, MSD, mucolipidosis, N-P (A/B), N-P (C), Sandhoff, SAS or TSD B1, comprising: determining a baseline level of a glucosidase in a tissue sample from the patient with the disorder before treatment with an agent; comparing a level of the a glucosidase in a tissue sample from the patient with the disorder after treatment with the agent with the baseline level; wherein a decrease relative to the baseline indicates a positive treatment outcome.
- 38. (Withdrawn) A method of monitoring a patient with Pompe's disease, comprising: determining a baseline level of a glucosidase in a tissue sample from the patient with the disorder before treatment with the agent; comparing a level of the a-glucosidase in a tissue sample from the patient after treatment with the agent with the baseline level; wherein an increase relative to the baseline indicates a positive treatment outcome.
- 39. (Previously Amended) A method of diagnosing or monitoring a lysosomal storage disorder in a patient subject, comprising:

obtaining a first sample from the patient subject;

measuring a first level of a saposin in the first sample obtained from the patient subject; comparing the first level to a baseline level, wherein the baseline level is the level of the saposin as determined in a control population of patients subjects unaffected by the lysosomal storage disorder;

determining a presence or extent of a lysosomal storage disorder when the first level is similar or different than the 95<sup>th</sup> percentile of the baseline level of at least the two saposins in the control population;

9

wherein,

## ATTORNEY DOCKET: MAYO-0005 (128675.00014)

- (i) the similarity of the first level compared to the baseline level is an indicator of absence of the lysosomal storage disorder in the patient subject;
- (ii) the difference of the first level compared to the baseline level is an indicator of presence or extent of the lysosomal storage disorder in the patient subject;
  - (iii) the saposin comprises saposin A, saposin B, saposin C, saposin D;
  - (iv) the first sample is plasma; and
  - (v) the baseline level and the first level are about equal to a percent elevation level for the lysosomal storage disorder <u>listed in Table 2</u>.

4643366v.2 10